**THUASNE®**

# Action Reliever®

<b>fr</b>	Genouillère souple de décharge .....	6
<b>en</b>	Soft offloading knee brace .....	8
<b>de</b>	Weiche Kniestütze .....	11
<b>nl</b>	Soepele ontlastende kniebrace .....	14
<b>it</b>	Ginocchiera di scarico morbida .....	16
<b>es</b>	Rodillera flexible de descarga .....	19
<b>pt</b>	Joelheira flexível de alívio do peso .....	22
<b>da</b>	Smidig, aflastende knæbind .....	25
<b>fi</b>	Kuormitusta keventävä taipuisa polvituki .....	27
<b>sv</b>	Flexibelt knästöd för avlastning .....	30
<b>el</b>	Μαλακή επιγονατίδα αποφόρτισης .....	32
<b>cs</b>	Odhlečovací chránič .....	35
<b>pl</b>	Elastyczna odciążająca orteza stawu kolanowego .....	38
<b>lv</b>	Elastīga atbrīvojoša ceļa ortoze .....	41
<b>lt</b>	Lankstus apkrovos mažinimo antkelis .....	43
<b>et</b>	Koormust vähendav elastne põlve ortoos .....	46
<b>sl</b>	Mehka razbremenilna kolenska opornica .....	48
<b>sk</b>	Odfahčovacia kolenná ortéza .....	51
<b>hu</b>	Puha tehermentesítő térdrögzítő .....	54
<b>bg</b>	Мека наколелка за разтоварване .....	56
<b>ro</b>	Genunchieră flexibilă pentru descărcare .....	59
<b>ru</b>	Мягкий наколенник для снятия напряжения .....	62
<b>hr</b>	Mekana ortoza, potpora za koljena .....	65
<b>zh</b>	减压用软质护膝 .....	67
<b>ar</b>	دعامة ركبة مرنة لتخفيف الثقل .....	71

Instructions d'utilisation - Instructions for use - Gebrauchsanweisung - Gebruiksaanwijzingen - Istruzioni per l'uso - Instrucciones de uso - Instruções de utilização - Anvisninger vedrørende brug - Käyttöohjeet - Råd vid användning - Οδηγίες χρήσης - Pokyny k používání - Szczegóły użytkowania - Lietošanas instrukcija - Naudojimo instrukcijos - Kasutusjuhised - Navodila za uporabo - Pokyny na používanie - Használati utasítás - Инструкции за употреба - Instrucțiuni de utilizare - Инструкция по эксплуатации - Upute za uporabu - 使用指南 - إرشادات الاستعمال

	(a)		(b)		
	cm	in	cm	in	
	1	38 - 41	15 - 16 ¼	32 - 34	12 ¼ - 13 ½
	2	42 - 45	16 ¼ - 17 ¾	35 - 37	13 ¾ - 14 ½
	3	45 - 48	17 ¾ - 19	38 - 41	15 - 16 ¼
	4	48 - 51	19 - 20 ½	42 - 45	16 ¼ - 17 ¾
	5	50 - 53	20 ½ - 21 ½	46 - 49	18 ¼ - 19 ¼
	6	53 - 56	21 - 22 ½	50 - 53	20 ½ - 21 ½
	7	56 - 61	22 - 24	53 - 56	21 - 22 ½
	8	61 - 66	24 - 26 ½	56 - 59	22 - 23 ¼

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## **SOFT OFFLOADING KNEE BRACE**

#### **Description/Destination**

The device is intended only for the treatment of the indications listed and for patients whose measurements correspond to the sizing table.

Please note that the orthosis can be used as follows:

- Medial arthrosis on the right/lateral arthrosis on the left,
- Medial arthrosis on the left/lateral arthrosis on the right.

Available in 8 sizes.

#### The device is composed of:

- an anatomically-shaped knitting integrating a light and soft knit at the back of the knee,
- 1 rigid upright and 1 flexible stay,
- 4 inelastic straps (❶❷❸❹),
- a patellar insert to support patella.

Depending on country of sale, additional accessories/spare parts could be available.

## **Composition**

Textile components: polyamide - elastane - polyester - silicone - polyurethane - polyethylene - cotton

Rigid components: steel - polyamide - polyurethane - silicone

## **Properties/Mode of action**

Adjustable 3-point leverage system with inelastic straps to provide alignment support and off-loading on the affected compartment.

An anatomically-shaped knitting for even compression.

The brace is maintained on the leg thanks to the silicone-coated threads at the top of the knee brace.

Easy to use thanks to the magnetics fasteners.

## **Indications**

Symptomatic unicompartimental femorotibial osteoarthritis (mild to moderate).

Knee off-loading for post-traumatic, post-operative or degenerative conditions.

Knee pain and/or swelling management.

## **Contraindications**

Do not use the product if the diagnosis has not been confirmed.

Do not use in the event of known allergy to any of the components.

Do not apply the product in direct contact with broken skin or an open wound without an adequate dressing.

Do not use in the event of major venous thromboembolic history without thromboprophylaxis.

Do not use if presence of implanted stimulation devices (as pacemaker,...).

## **Precautions**

Verify the product's integrity before every use.

Do not use the device if it is damaged.

Choose the appropriate size to fit the patient, referring to the size chart.

If the measurements taken do not correspond to the same size, favor the size corresponding to measurement **(a)**.

It is recommended that a healthcare professional supervises the first application.

Strictly comply with your healthcare professional's prescription and recommendations for use.

This product is intended for the treatment of a given condition. Its duration of use is limited to this treatment only.

For hygiene and performance reasons, do not re-use the product for another patient.

It is recommended to adequately tighten the device to achieve support/immobilisation without restricting blood circulation.

In case of history of venous or lymphatic disorders, consult a healthcare professional.

In the event of discomfort, significant hindrance, pain, variation in limb volume, abnormal sensations or change in colour of the extremities, remove the device and consult a healthcare professional.

In the event of any modification in the product's performance, remove it and consult a healthcare professional.

Before any sports activity, check the compatibility of the use of this medical device with your healthcare professional.

This device has magnetic components.

Do not wear the device in a medical imaging machine.

Do not use the device in case of application of certain products on the skin (creams, ointments, oils, gels, patches...).

## **Undesirable side-effect**

This device can cause skin reactions (redness, itching, burns, blisters, etc.) or wounds of various degrees of severity.

Any serious incidents occurring related to the device should be reported to the manufacturer and to the competent authority of the Member State in which the user and/or patient is resident.

## Instructions for use/Application

It is recommended to wear the device directly on the skin, unless contraindicated. Before you put it on, first loosen all the straps.

### Important:

- When fitting the orthosis, always make sure that the two intersecting self-gripping straps (crossover self-gripping straps) ❶ are always positioned centrally during adjustment: check the spacing and position of the two crossover self-gripping straps.

The length of the two ends of each crossover self-gripping strap protruding from the crossover strap cover must be identical.

- To avoid damage to the knit fabric, please make sure that the self-gripping pads are always fastened to the designated self-gripping surfaces when putting it on and taking it off.

### First fitting (straps adjustment):

Have the patient sit on the edge of a chair, and have the patient bend his/her knee to a position of approximately 30-60° of flexion.

- Slip the knee support on the leg with your foot extended.  
To make it easier to put on, use the pull tabs at the top of the knee brace.
- Make sure that your knee cap is positioned in the centre of the patellar insert.

### **Start the adjustment of the two rear straps. First fit the lower leg ❷ and then the upper leg ❸ strap ❶.**

- To adjust the length of the two rear straps, it is necessary to adjust the length of the end opposite the self-gripping strap intersection which has a straight cut edge ❷.
- To adapt the strap, remove the self-gripping pad and shorten the strap to the desired length ❸.  
Then re-attach the self-gripping pad and fix it on the strap.
- The self-gripping pad should be positioned centrally on the length of the strap so that the tension can be increased if necessary.
- If creases form during strap adjustment, the angle of each strap can be adjusted and positioned individually by opening the two rear hook and loop openings on the crossover self-gripping strap (round end) ❹.

### **Then adjust the two front straps. First fit the lower leg ❹ and then the upper leg ❺ strap ❶.**

- To adjust the length of the two front straps, it is necessary to adapt the length of the end of the strap which is directly located on the crossover self-gripping strap ❶.
- To adapt the strap, remove the self-gripping pad and shorten the strap to the desired length ❺.

Take care to not cut any strap too short.

- The self-gripping pad should preferably be positioned centrally on the length or near the fastening loop with the magnetic fastener so that the tension can be increased if necessary.
- To close the strap, grasp into the loop at the end of the strap ❻ and pull the end of the strap towards the self-gripping pad with the magnetic fastener so that the magnetic fastener can find itself and close automatically.
- If creases form during strap adjustment, the angle of each strap can be adjusted and positioned individually by opening the two rear hook and loop openings on the crossover self-gripping strap (round end) ❽.
- To change the tension while wearing the orthosis, it is only necessary to change the position of the self-gripping ends of the straps on the straps.
- If the lower end of the device presses into the calf, the hinged upright can be bent outwards at the distal end appropriately.

### Before you take it off:

To loosen the strap system, please open the magnetic fasteners of the two front straps.

It is not necessary to readjust the rear straps ❷ and ❸.

Then grasp into the orthosis from below and pull it off downwards.

## Care/Maintenance

Close the self-fastening tabs before washing.

Remove only the rigid upright before washing. To do this:

- Pull on the top of the product to disengage the upper end of the lateral uprights from of its casing.
- Take note of its original position (flexed towards the back), in order to reposition it correctly.

It is recommended to remove the patellar insert before washing. Take note of its original position (flat area outside, wide area above), in order to reposition it correctly. Machine washable at 30°C (delicate programme). If possible use a washing net. Do not use detergents, fabric softeners or aggressive products (products containing chlorine). Squeeze out excess water. Dry away from any direct heat source (radiator, sun, etc.). If the device is exposed to seawater or chlorinated water, make sure to rinse it in clear water and dry it.

## Storage

Store at room temperature, preferably in the original packaging.

## Disposal

Dispose of in accordance with local regulations.

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### **Justera därefter de två främre remmarna. Sätt först fast den nedre remmen © och därefter den övre remmen ①.**
**THUASNE UK Ltd**
Thuasne SAS - SIREN/RCS Nanterre 542 091 186