


THUASNE

BIPLAST

fr	Bande adhésive de contention	04
en	Adhesive compression bandage	07
de	Haftender Kompressionsverband	09
nl	Zelfklevend steunverband	12
it	Fascia adesiva di contenimento	15
es	Venda adhesiva de contención	18
pt	Banda adesiva de contenção	20
fi	Liimaside	23
sv	Självhäftande stödband	26
el	Αυτοκόλλητη ταινία συγκράτησης	28
pl	Pas podtrzymujący samoprzylepny	31
lv	Līpoša balsta saite	34
sk	Adhezívny fixačný obvaz	37
hu	Öntapadós rugalmas fásli	39
bg	Адхезивен бинт за задържане	42
ro	Bandă adezivă de compresie	45
ru	Бинт липкий фиксирующий	48
ko	접착식 고정 붕대	51
ar	رباط ضاغط لاصق	53

MD

fr	Dispositif médical
en	Medical device
de	Medizinprodukt
nl	Medisch hulpmiddel
it	Dispositivo medico
es	Producto sanitario
pt	Dispositivo médico
fi	Lääkinnällinen laite
sv	Medicinteknisk produkt
el	Ιατροτεχνολογικό προϊόν
pl	Wyrób medyczny
lv	Medicīniska ierīce
sk	Zdravotnícka pomôcka
hu	Orvostechikai eszköz
bg	Медицинско изделие
ro	Dispozitiv medical
ru	Медицинское изделие
ko	의료 기기
ar	جهاز طبي

A

fr **BANDE ADHÉSIVE DE CONTENTION**

Description

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ADHESIVE COMPRESSION BANDAGE

Description

Short-stretch adhesive compression bandage.

Composition:

Viscose, polyamide, paper, rubber.

Properties/Mechanism of action

Adhesive compression bandage.

Aerated elastic fabric.

Multidirectional chevron patterned adhesive backing.

Short stretch < 70%.

Protective, release backing paper.

Biplast bandages are short-stretch bandages indicated in the treatment of trauma, vascular and lymphatic conditions.

Use:

Biplast is intended to be placed on the body part or a joint affected by the patient's condition and in accordance with one of the intended medical indications.

Indications

Soft tissue or musculoskeletal injuries.

Taping or strapping.

Fracture of small bone (finger, foot, etc.).

Bone deformity.

Vascular or lymphatic conditions requiring short-stretch or multicomponent bandaging.

Biplast bandages can be used to secure dressings and splints.

Instructions for use/Application

Follow the instructions given to you by a health care professional on the fitting method and on how long the bandage is to be worn.

Select the desired bandage length.

Remove the protective backing paper **A**.

Contraindications

Do not use Biplast bandages in the following cases:

Do not apply the product in direct contact with broken skin or an open wound without an adequate dressing.

Do not use in the event of known allergy to any of the components.

Do not use Biplast bandages in the following cases if they are being used for a vascular or lymphatic condition:

Phlegmasia cerulea dolens (painful blue edema with arterial compression).

Septic thrombosis.

Do not use Biplast bandages in the following cases if they are being used for a vascular or lymphatic condition of the lower limbs:

Peripheral arterial disease (PAD) with an ABI < 0.6.

Advanced diabetic microangiopathy.

Decompensated heart failure.

Precautions

Ask your doctor for advice if your skin is particularly sensitive.

In the event of discomfort, significant hindrance, pain, variation in limb volume, abnormal sensations or change in color of the extremities, remove the device and consult a health care professional.

Follow the instructions given to you by a health care professional on

fitting and on the recommended treatment duration.

Do not wear the bandage for an extended period of time without getting a medical opinion.

Make sure you cover any infected or open skin wound before using the bandage.

Any serious incidents occurring related to the device should be reported to the manufacturer and to the competent authority of the Member State in which the user and/or patient is resident.

Do not re-use used bandages.

Adverse effects

This device may cause skin reactions such as pressure marks, redness, itching, burns, and even blisters.

Care/Maintenance

Single-use, disposable device.

Store the device at room temperature in its original packaging.

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©Thuasne - 2019601 (2020-10)